

BIPOLAR CABLES



Timesco Bipolar Cables are supplied with 2 universal "Banana Type" plugs for fitting to most popular electrosurgery machines. Each cable also has a tag to allow the user to easily record the number of uses, so that the cable may be replaced within its safe working life.

Timesco Bipolar Cables are made from supple elastic silicon which will not crinkle or crack under normal use.



Two types of cable are available both 3 metres in length: The European fitting "square end" and the USA fitting "round end". (See forceps for these types illustrated on pages 2 and 3 of this publication.) Standard machine ends are supplied on both cables.

Timesco Bipolar Forceps and Cables are fully steam autoclaveable. Ultrasonic cleaning with volatile vapourising solvents is not permissible. Care and maintenance instructions are available with each cable.

In this brochure only a selection of the popular Bipolar products are illustrated. Our range is very extensive and if the items required are not illustrated please request further details

87.875.20
TIMESCO BIPOLAR CABLE
3 METRES USA SILICON



87.850.20
TIMESCO BIPOLAR CABLE
3 METRES EURO SILICON



BIPOLAR FORCEPS AND ACCESSORIES

For samples or further information on any product, please contact our sales team on (020) 7511 1234.



Timesco
SURGICAL + MEDICAL
INSTRUMENTS FOR LIFE



Timesco of London, Timesco House, 1 Knights Road, London E16 2AT.
Telephone: +44 (020) 7511 1234 Fax: +44 (020) 7511 7888 Email: surgical@timesco.com
www.timesco.com

E & O.E. Copyright October 2002. Ref B.S.001/1

'Timesco' and 'Instruments for Life' are registered trademarks of Timesco of London. All rights reserved.

Timesco
SURGICAL + MEDICAL
INSTRUMENTS FOR LIFE

BIPOLAR FORCEPS - EUROPEAN STYLE



We produce our forceps with the popular "European Fitting" for connecting to bipolar cables.

Timesco Bipolar forceps and accessories are precision made instruments which are specially coated with a non-cracking elastic insulation coating which is guaranteed to be 100% safe to use for the operator and the patient.

Timesco Bipolar Forceps are made from durable stainless steels and industry standard insulation coatings.

The forceps will perform equally well in wet or dry operative fields. The coating used on our forceps is double thickness which is leak tested electronically to at least 40,000 volts. We are one of the few companies who can offer this guarantee.

STANDARD FORCEPS STRAIGHT



87.640.20 STRAIGHT 6" (15CMS)
FINE TIP 0.50MM

ALSO AVAILABLE:
87.645.20 STRAIGHT 6" (15CMS) BROAD TIP 1.0MM
87.655.20 STRAIGHT 7" (17.5CMS) FINE TIP 0.50MM
87.660.20 STRAIGHT 7" (17.5CMS) BROAD TIP 1.0MM
87.665.20 STRAIGHT 7.75" (19.5CMS) FINE TIP 0.50MM
87.668.20 STRAIGHT 7.75" (19.5CMS) BROAD TIP 1.0MM
87.670.20 STRAIGHT 8.75" (21.5CMS) FINE TIP 0.50MM
87.672.20 STRAIGHT 8.75" (21.5CMS) BROAD TIP 1.0MM


STANDARD FORCEPS ANGLED



87.640.40 ANGLED 6" (15CMS)
FINE TIP 0.50MM

ALSO AVAILABLE:
87.645.40 ANGLED 6" (15CMS) BROAD TIP 1.0MM
87.655.40 ANGLED 7" (17.5CMS) FINE TIP 0.50MM
87.660.40 ANGLED 7" (17.5CMS) BROAD TIP 1.0MM
87.665.40 ANGLED 7.75" (19.5CMS) FINE TIP 0.50MM
87.668.40 ANGLED 7.75" (19.5CMS) BROAD TIP 1.0MM
87.670.40 ANGLED 8.75" (21.5CMS) FINE TIP 0.50MM
87.672.40 ANGLED 8.75" (21.5CMS) BROAD TIP 1.0MM

87.570.20 JEWELLER
BIPOLAR FORCEPS STRAIGHT 4.5" (11.25CMS)
MICRO TIP 0.25MM



87.570.40 JEWELLER
BIPOLAR FORCEPS CURVED 4.5" (11.25CMS)
MICRO TIP 0.25MM



87.580.20 MCPHERSON
BIPOLAR FORCEPS STRAIGHT 4" (10CMS)
FINE TIP 0.50MM




87.580.40 MCPHERSON
BIPOLAR FORCEPS ANGLED 4" (10CMS)
FINE TIP 0.50MM




BIPOLAR BAYONET FORCEPS

ALSO AVAILABLE:
87.678.40 BAYONET 7" (17.5CMS) BROAD TIP 1.0MM
87.680.20 BAYONET 7.75" (19.5CMS) FINE TIP 0.50MM
87.680.40 BAYONET 7.75" (19.5CMS) BROAD TIP 1.0MM



87.678.20 BAYONET 7" (17.5CMS)
FINE TIP 0.50MM

ALSO AVAILABLE:
87.720.20 ANGLED UP 7" (17.5CMS) FINE TIP 0.50MM
87.720.40 ANGLED UP 7" (17.5CMS) BROAD TIP 1.0MM
87.722.20 ANGLED UP 7.75" (19.5CMS) FINE TIP 0.5MM
87.722.40 ANGLED UP 7.75" (19.5CMS) BROAD TIP 1.0MM



87.730.20 ANGLED DOWN 7" (17.5CMS)
FINE TIP 0.5MM

ALSO AVAILABLE:
87.730.40 ANGLED DOWN 7" (17.5CMS) BROAD TIP 1.0MM
87.732.20 ANGLED DOWN 7.75" (19.5CMS) FINE TIP 0.5MM
87.732.40 ANGLED DOWN 7.75" (19.5CMS) BROAD TIP 1.0MM


BIPOLAR FORCEPS - USA STYLE

We also produce Bipolar Forceps with the USA style "Two Pin" fitting. This is more widely used not only in the USA but also in South America and the Far East. A selection of these are shown in this brochure. Our range of these patterns is also quite extensive, if further details are required please contact us.



USA style
"Two Pin" fitting.


87.865.20 JEWELLER
FORCEPS STRAIGHT 4.5" (11.25CMS)
TIP 0.25MM USA




87.855.20 ADSON FORCEPS
STRAIGHT 4.75" (11.75CMS)
TIP 1.0MM USA



87.860.20 CUSHING FORCEPS
STRAIGHT 6.75" (16.75CMS)
TIP 1.0MM USA



87.862.20 GERALD FORCEPS
BAYONET 7.25" (18CMS)
TIP 1.0MM USA



87.868.20 SEMKIN FORCEPS
STRAIGHT 5.5" (13.75CMS)
TIP 1.0MM USA

Timesco Bipolar Forceps are coated with a special Blue coating to help distinguish them from "blackened" or "laser" instruments. This helps avoid confusion in the operating theatre as some older versions of electrosurgery instruments are still available with a black coating.