


SCISSORS, PLASTER SAWS & KNIVES

We have also introduced a special selection of instruments with Tungsten Carbide jaws. Due to their incredibly hard properties, Tungsten Carbide jaws will not wear down or blunt as quickly as normal steel jaws, even during repeated use. Tungsten Carbide jaws also benefit from being very sharp without becoming too brittle. We have also added a serrated edge to a selection of our Tungsten Carbide instruments to prevent slippage, especially useful when cutting synthetic materials.


43.025.50
TIMESCO BANDAGE SCISSORS 7" (17.5CMS)
WITH EXCENTRIC JOINT, ONE SERRATED BLADE



12.485.15
LISTER BANDAGE SCISSORS 7.5" (18.75CMS)
WITH PROBE POINT



12.955.20
LISTER BANDAGE SCISSORS 7" (17.5CMS)
WITH TUNGSTEN CARBIDE JAWS



12.993.10
TIMESCO UNIVERSAL SCISSORS
5" (12.5CMS)
WITH TUNGSTEN CARBIDE JAWS
AND UPPER SERRATED JAW
FOR SYNTHETIC MATERIALS



Our new Universal Cutter with Tungsten Carbide Jaws also has a specially serrated blade to prevent slippage and aid the cutting of synthetic materials.



43.110.10
ENGEL PLASTER SAW 6" (15CMS)

43.110.20
SPARE BLADE ONLY
FOR ENGEL SAW



12.695.25
UNIVERSAL SCISSORS 4.75" (11.75CMS)
WITH SERRATED BLADE

43.040.25
ESMARCH PLASTER KNIFE 7" (18CMS)



For samples or further information on any product, please contact our sales team on (020) 7511 1234.



Timesco
SURGICAL + MEDICAL
INSTRUMENTS FOR LIFE



Timesco of London, Timesco House, 1 Knights Road, London E16 2AT.
Telephone: +44 (020) 7511 1234 Fax: +44 (020) 7511 7888 Email: surgical@timesco.com
www.timesco.com

E & O.E. Copyright October 2002. Ref B.S.002/1

'Timesco' and 'Instruments for Life' are registered trademarks of Timesco of London. All rights reserved.



PLASTER INSTRUMENTS

Timesco
SURGICAL + MEDICAL
INSTRUMENTS FOR LIFE

PLASTER SHEARS

We produce one of the finest ranges of quality surgical instruments in the World. Our extensive Plaster Instruments range is no exception. Each item is made to offer the user an instrument which has been developed to provide the highest level of performance and durability. The nature of Plaster Instruments requires them to be robust enough to withstand repeated heavy use. As with all Timesco reusable instruments, our Plaster Instruments are manufactured from the highest quality surgical stainless steel, are fully autoclaveable, and Guaranteed for Life against faulty manufacture and materials.

12.710.05
TIMESCO UTILITY SHEAR 7.5" (18.75CMS)
WITH GREEN PLASTIC HANDLES



ALSO AVAILABLE:
12.710.20 BLUE
12.710.30 BLACK
12.710.35 RED

43.130.10
BOHLER (BRUN) PLASTER SHEAR 10" (25CMS)
WITH ONE SERRATED BLADE



43.130.20
BOHLER PLASTER SHEAR 10" (25CMS)
WITH SMOOTH BLADES




43.190.20
TIMESCO PLASTER SHEAR
9.5" (23.75CMS)
WITH TUNGSTEN CARBIDE JAWS
AND UPPER SERRATED JAW
FOR SYNTHETIC MATERIALS



43.140.10
GUY PLASTER SHEAR 11.5" (29CMS)



43.160.10
SEUTIN PLASTER SHEAR 9" (23CMS)



ALSO AVAILABLE
43.160.20
SEUTIN PLASTER SHEAR 8" (20CMS)
WITH LOCK

PLASTER SHEARS, SPREADERS AND BREAKERS

ALSO AVAILABLE:
43.180.20
STILLE PLASTER SHEAR 13" (32CMS)
WITH SERRATED UPPER BLADE
43.180.30
STILLE PLASTER SHEAR 15" (38CMS)
WITH SERRATED UPPER BLADE




43.180.10
STILLE PLASTER SHEAR 9" (23CMS)
WITH SERRATED UPPER BLADE



43.250.10
TRI-BLADE PLASTER SPREADER
DOUBLE SPRING 9" (23CMS)

43.270.20
WOLFF PLASTER BREAKER 9.5" (24CMS)



ALSO AVAILABLE:
43.270.10
WOLFF PLASTER BREAKER 7" (18CMS)



43.230.10
HENNIG PLASTER SHREDDER 11" (28CMS)



43.210.10
DAWS PLASTER SPREADER 18" (45CMS)



The jaw of the Tri Blade spreader has been specially designed to be fine enough to be inserted between plaster and synthetic casts which have been cut with an electric saw blade.



Our new Synthetic Shears with Tungsten Carbide Jaws also have a specially serrated blade to prevent slippage and aid the cutting of synthetic materials.